## **Specimen Collection and Transportation Information**

Test Name	Mnemonic	Section	<b>Specimen Type</b>
Urine Sugar Chromatography	SUCHR	Clinical Chemistry	Urine

- Preparation of the Patient: Not required
- Type of Collection Container: Sterile urine container without any preservative
- Types and Amounts of Preservatives or Anticoagulants: Not required
- Quantity of Specimen to be Collected: 2-5 ml of urine
  - o Minimum Quantity: 1.0 ml (1000 μl)
- Transportation Temperature: Not Required
- Need for Special Timing for Collection (where applicable): Ideally urine should be collected at a time when symptoms are most severe
- Specimen Stability Information: As same as for the routine urine specimen
- **Special Instruction:** Early morning samples preferred, otherwise random urine in a container with no preservatives
- Rejection Criteria: Quantity Not Sufficient (QNS)
- Need for Appropriate Clinical Data, When Indicated (Patient history): Not Required

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Thin layer chromatography	After 2 days	on Tuesday	12:01 AM	12:01 AM	12:01 AM