Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Amphetamine	AMPH	Clinical Chemistry	Urine

- Preparation of the Patient: Freshly Random Urine Sample Required
- Type of Collection Container: Urine Container
- Types and Amounts of Preservatives or Anticoagulants: NA
- Quantity of Specimen to be Collected: 20 ml
 - o **Minimum Quantity:** 4-5ml
- Transportation Temperature: 2-8°C After 48 Hours 20°C
- Need for Special Timing for Collection (where applicable): NA
- Specimen Stability Information: NA
- Special Instruction: NA
- Rejection Criteria: Wrong Container, QNS, Wrong Identification
- Need for Appropriate Clinical Data, When Indicated (Patient history): Not Required

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Enzyme Multiplied Immunoassay Technique (EMIT)	Same Day	Working Days	11:59 pm	12:00 am	12:00 am