



Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Barbiturates	BARBI	Clinical Chemistry	Urine

- **Preparation of the Patient:** Freshly Random Urine Sample Required
- **Type of Collection Container:** Urine Container
- **Types and Amounts of Preservatives or Anticoagulants:** NA
- **Quantity of Specimen to be Collected:** 20 ml
 - **Minimum Quantity:** 4-5ml
- **Transportation Temperature:** 2-8°C After 48 Hours - 20°C
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** NA
- **Special Instruction:** NA
- **Rejection Criteria:** Wrong Container, QNS, Wrong Identification
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Enzyme Multiplied Immunoassay Technique (EMIT)	Same Day	Working Days	11:59 pm	12:00 am	12:00 am