



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Breath Hydrogen	BHT	Clinical Chemistry	Breath

- **Preparation of the Patient:** Patient should be in fasting
- **Type of Collection Container :** Syringe
- **Types and Amounts of Preservatives or Anticoagulants :** NA
- **Quantity of Specimen to be Collected:** Breath 30cc
 - **Minimum Quantity:** 15cc
- **Transportation Temperature :** Room Temperature
- **Need for Special Timing for Collection (where applicable):** Appointment Required
- **Specimen Stability Information:** NA
- **Special Instruction:** NA
- **Rejection Criteria:** NA
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Gas Chromatography	Same Day	Working Days	11:00 am	NA	NA