



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital

Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Lithium	LIT	Clinical Chemistry	Serum

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Plain Tube or Gel tube
- **Types and Amounts of Preservatives or Anticoagulants:** NA
- **Quantity of Specimen to be Collected:** Serum Or Clotted Sample 3-5cc
 - **Minimum Quantity:** 3cc
- **Transportation Temperature:** 2-8°C After 48 Hours - 20°C
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** NA
- **Special Instruction:** NA
- **Rejection Criteria:** Wrong Tube, QNS, Hemolysis, Wrong Identification
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Ion Selective Electrode (ISE) method	Same Day	Working Days	12:00 pm	12:00 am	12:00 am