Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Porphyrins	PORPHYRINS	Clinical Chemistry	Urine

• Preparation of the Patient: NA

• Type of Collection Container: Freshly urine sample required

Types and Amounts of Preservatives or Anticoagulants: NA

• Quantity of Specimen to be collected: 50ml

o Minimum Quantity: 20 ml

• Transportation Temperature: 18-25°C (For 8 Hour) Or Keep at 2-8°C and used within 24 hours

• Need for Special Timing for Collection (where applicable): NA

• Specimen Stability Information: NA

• **Special Instruction:** Available

• Rejection Criteria: Wrong Container, QNS, Wrong Identification

• Need for Appropriate Clinical Data, When Indicated (Patient history): Not Required

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Chemical Method	Next Day	Every day	8:00 am	8:00 am	8:00 am