



Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Sweat Chloride	SWCL	Clinical Chemistry	Sweat

- **Preparation of the Patient:** Patient well hydrated, free of acute illness, and not receiving mineral corticoids.
- **Type of Collection Container:** Macro duct
- **Types and Amounts of Preservatives or Anticoagulants :** NA
- **Quantity of Specimen to be Collected:** Sweat Sample 6-10ul
 - **Minimum Quantity:** 6ul
- **Transportation Temperature:** Liquid sweat from macro duct can be stored in sealed macro duct tubing for up to 72 hours at 4°C
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** NA
- **Special Instruction:** NA
- **Rejection Criteria:** Minimum sample volume required 6 to 10 micro liters. Less than stated volume test could not be performed
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
conductivity	Same Day	Working Days	2:00 pm	NA	NA