



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Acute Leukemia Comprehensive/Basic Panel or Lymphoma Basic Panel	ALCP/ALBP OR LBP	Histopathology	Blood / Bone Marrow/ Body Fluid

- **Preparation of the Patient:** N/A
- **Type of Collection Container:** EDTA and Sodium Heparin Tube (purple and green top)
- **Types and Amounts of Preservatives or Anticoagulants:** EDTA/Sodium heparin
- **Quantity of Specimen to be Collected:** 4ml
 - **Minimum Quantity:** 3ml
- **Transportation Temperature:** 2-8°C
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** Temperature should be maintained
- **Special Instruction:** As the sample is collected it should be transported as soon as possible
- **Rejection Criteria:** QNS, Hemolysed
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** % of blast cells, lymphocytosis? TLC count, CBC report

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
immunophenotyping	3 working days	Mon- Fri	3:00pm	3:00pm	3:00pm