Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type	
Amoebiasis	AMOEBA	Clinical Microbiology	Serum or Clotted Blood	

• Preparation of the Patient: NA

• Type of Collection Container: Gel tube

Types and Amounts of Preservatives or Anticoagulants: NA

• Quantity of Specimen to be Collected: 3-5cc

o Minimum Quantity: 3cc

- Transportation Temperature: If delay in between 2-4hrs at Room temp, If delay more than 2 hrs at refrigerated temp
- Need for Special Timing for Collection (where applicable): NA
- Specimen Stability Information: If delay in between 2-4hrs at Room temp, If delay more than 2 hrs at refrigerated temp

• Special Instruction: NA

• Rejection Criteria: Hemolysed blood not accepted

• Need for Appropriate Clinical Data, When Indicated (Patient history): Not required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
ELISA	Wednesday & Friday	Tuesday & Thursday	00.01	10.00	Contact nearest Collection Point