Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Anti-Streptolysin-O Titre	ASO	Clinical Microbiology	Serum

• Preparation of the Patient: NA

• Type of Collection Container: Gel tube

• Types and Amounts of Preservatives or Anticoagulants: NA

• Quantity of Specimen to be Collected: 3-5cc

Minimum Quantity: 3cc

• Transportation Temperature: Room temp

• Need for Special Timing for Collection (where applicable): NA

• Specimen Stability Information: If delay in between 2-4hrs at Room temp, If delay more than 4 hrs. at refrigerated temp

• Special Instruction: NA

• Rejection Criteria: Hemolyzed blood not accepted

• Need for Appropriate Clinical Data, When Indicated (Patient history): Not required

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Agglutination	MON thru FRI, same day	Daily	00.01	10.00	Contact nearest Collection Point