



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital  
Clinical Laboratories

### Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Bone	BO	Histopathology	Tissue Bone Biopsy Sample

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Wide mouth Plastic Container
- **Types and Amounts of Preservatives or Anticoagulants:** 10% buffered formalin 10 times volume of the specimen
- **Quantity of Specimen to be Collected:**
  - **Minimum Quantity:** As per physician & surgeon of the patient
- **Transportation Temperature:** Room Temperature
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** Should be immerse in formalin
- **Special Instruction:** Radiological Films
- **Rejection Criteria:** Without Radiological Films/Without formalin/Without history
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Complete clinical history

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Histology	10 working days	as per cut off	03:00 pm	NA	Contact Nearest Collection Point