Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type	
Enzyme Histochemistry	EHC	Histopathology	Fresh Muscle Biopsy Specimen	

- Preparation of the Patient: NA
- Type of Collection Container: Wide mouth Plastic Container
- Types and Amounts of Preservatives or Anticoagulants: NA
- Quantity of Specimen to be Collected:
 - o Minimum Quantity: As per surgeon
- Transportation Temperature: Room Temperature
- Need for Special Timing for Collection (where applicable): NA
- Specimen Stability Information: Send Immediately to the lab
- **Special Instruction:** Unfix, Fresh specimen
- Rejection Criteria: Received with fixative
- Need for Appropriate Clinical Data, When Indicated (Patient history): Complete clinical history

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Cryosectioning	15 working days	As per cut off	03:00 pm	NA	NA