Specimen Collection and Transportation Information

Test Name	Test Name Mnemonic		Specimen Type	
Fine Needle Aspiration Cytology	FNACP	Histopathology	Aspirated Fluid Samples / Slides	

- Preparation of the Patient: Prepare patient for procedure as per policy
- Type of Collection Container: Slides in slide mailer and sample in syringe
- Types and Amounts of Preservatives or Anticoagulants: Aspiration in alcohol and spray (95% alcohol) fixed slides
- Quantity of Specimen to be Collected: Aspirate as material as possible
 - o Minimum Quantity: NA
- Transportation Temperature: Room temperature
- Need for Special Timing for Collection (where applicable): NA
- Specimen Stability Information: NA
- **Special Instruction:** Slides should be properly / adequately prepared (3 alcohol fixed and 1 dry) for each attempt
- Rejection Criteria: Unsatisfactory Smear preparation or insufficient quantity
- Need for Appropriate Clinical Data, When Indicated (Patient history): Complete clinical information, previous surgeries etc.

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Cytology	after 04 working days	all Days except Sunday	Daily before 4:00 pm and Public holidays (before 3:00 pm)	NA	Coordinate with stat labs