



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Gross Only	GO	Histopathology	Tissue Biopsy Sample or Other Nonspecific Sample

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Wide mouth Plastic Container
- **Types and Amounts of Preservatives or Anticoagulants:** 10% buffered formalin 10 time volume of specimen (in case of tissues In container)
- **Quantity of Specimen to be Collected:**
 - **Minimum Quantity:** As per surgeon
- **Transportation Temperature:** Room temperature
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** NA
- **Special Instruction:** Please see the specimen Gross only list
- **Rejection Criteria:** NA
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Complete clinical history along with sample

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
NA	same day	as per cut off	12:00 pm	NA	Contact Nearest Collection Point



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Gross Only	GO	Histopathology	Tissue Biopsy Sample or Other Nonspecific Sample