



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital  
Clinical Laboratories

### Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
HLADRB1 BY PCR	HLA DR-2 BY PCR	Molecular Pathology	Whole blood in EDTA

- **Preparation of the Patient:** NA
- **Type of Collection Container:** 04 ml EDTA Tubes(Purple Top)
- **Types and Amounts of Preservatives or Anticoagulants:** EDTA (7.2mg)
- **Quantity of Specimen to be Collected:** 10ml Blood in EDTA
  - **Minimum Quantity:** NA
- **Transportation Temperature:** 2-8° C
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** 2-8°C
- **Special Instruction:** NA
- **Rejection Criteria:**
  - a) Hemolysis/ clotted
  - b) Leaky Container
  - c) In appropriate container
  - d) Improper labeled specimen
  - e) In adequate volume of specimen
  - f) Wrong source
  - g) Improper transportation temperature
  - h) Contaminated specimen with biological hazardous material
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** NA  
Patient detail required (HLA report copy or MR of L#)

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
PCR	Friday of every week	Friday of every week	0.01	10:00	Contact Nearest Collection Points



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital  
Clinical Laboratories

**Specimen Collection and Transportation Information**

Test Name	Mnemonic	Section	Specimen Type
HLADRB1 BY PCR	HLA DR-2 BY PCR	Molecular Pathology	Whole blood in EDTA