



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital  
Clinical Laboratories

**Specimen Collection and Transportation Information**

Test Name	Mnemonic	Section	Specimen Type
Gross Micro Large Specimen	LB	Histopathology	Tissue Biopsy Sample

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Wide mouth Plastic Container
- **Types and Amounts of Preservatives or Anticoagulants:** 10% buffered formalin 10 time volume of specimen
- **Quantity of Specimen to be Collected:** As per physician & surgeons
  - **Minimum Quantity:** As per physician & surgeons
- **Transportation Temperature:** Room temperature
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** Should be immerse in formalin
- **Special Instruction:** NA
- **Rejection Criteria:** Without formalin/Without history
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Complete clinical history along with sample

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Histology	06 working days	as per cut off	03:00 pm	NA	Contact Nearest Collection Point