Specimen Collection and Transportation Information

Test Name Mnemonic		Section	Specimen Type	
Latex Particle Agglutination	LPA	Clinical Microbiology	CSF (Cerebrospinal Fluid)	

• Preparation of the Patient: NA

• Type of Collection Container: Sterile container

Types and Amounts of Preservatives or Anticoagulants: NA

• Quantity of Specimen to be Collected: 1-3ml

Minimum Quantity: 0.5ml

• Transportation Temperature: Room temperature

• Need for Special Timing for Collection (where applicable): NA

• **Specimen Stability Information:** If delay in between 2-4hrs at Room temp, If delay more than 4 hrs. at refrigerated temp

• Special Instruction: NA

• Rejection Criteria: Unsterile container / sample in EDTA

• Need for Appropriate Clinical Data, When Indicated (Patient history): Required

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Agglutination	Mon to Fri same day	Daily, same day (Within 2hrs.)	00.01	10.00	Contact nearest Collection Point