**Specimen Collection and Transportation Information** 

Test Name Mnemonic		Section	Specimen Type	
Leptospira Antibodies	LAT	Clinical Microbiology	Serum / Clotted Blood	

• Preparation of the Patient: NA

• Type of Collection Container: Gel tube

Types and Amounts of Preservatives or Anticoagulants: NA

• Quantity of Specimen to be Collected: 3-5cc

o Minimum Quantity: 3cc

• Transportation Temperature: Room temperature

Need for Special Timing for Collection (where applicable): NA

• **Specimen Stability** Information: If delay in between 2-4hrs at Room temp, If delay more than 4 hrs. at refrigerated temp

• Special Instruction: NA

• Rejection Criteria: Hemolysed blood not accepted

• Need for Appropriate Clinical Data, When Indicated (Patient history): Not Required

Methodology	Reporting	Day	Cut Off Time	Cut Off Time	Cut Off Time Out
	Scheme	Performed	Main Lab	Karachi Points	of Karachi Points
Agglutination	Same day if received before 12:00 noon	Daily	00.01	10.00	Contact nearest Collection Point