



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Alcohol	ALC	Clinical Chemistry	Plasma

- **Preparation of the Patient:** NA
- **Type of Collection Container:** EDTA
- **Types and Amounts of Preservatives or Anticoagulants:** K2 EDTA
- **Quantity of Specimen to be collected:** 3cc blood in EDTA, collect without spirit swab
 - **Minimum Quantity:** 3cc
- **Transportation Temperature:** 2-8°C After 48 Hours – 20°C
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** NA
- **Special Instruction:** NA
- **Rejection Criteria:** : Wrong Tube, QNS, Hemolysis, Wrong Identification
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Enzymatic Reaction	Next Day	Daily	11:59 PM	12:00am	12:00am