



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Acetyl Choline Receptor Ab	ACETYL	Clinical Chemistry	Serum

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Plain tube or Gel Tube (yellow top)
- **Types and Amounts of Preservatives or Anticoagulants:** NA
- **Quantity of Specimen to be Collected:** 3-5cc
 - **Minimum Quantity:** 3cc
- **Transportation Temperature:** 18-25°C for 8 Hours, 2-8°C after 8 hours to 14 days
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** Assay within 8 hours store sample at room temperature (18-25°C),for 14 days store at 2-8°C
- **Special Instruction:** NA
- **Rejection Criteria:** Wrong Tube, QNS, Hemolysis, Wrong identification, Lipemic.
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Enzyme Linked Immunosorbent Assay (ELISA)	After 3 days	Friday	12:01 am	12:01 am	12:01 am



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Acetyl Choline Receptor Ab	ACETYL	Clinical Chemistry	Serum