



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Cyclosporin	Cyclo	Clinical Chemistry	Whole Blood

- **Preparation of the Patient:** NA
- **Type of Collection Container:** EDTA
- **Types and Amounts of Preservatives or Anticoagulants:** K2 EDTA
- **Quantity of Specimen to be Collected:** 2-3cc EDTA blood (whole blood)
 - **Minimum Quantity:** 3cc
- **Transportation Temperature:** 18-25°C (For 8 Hour) 2-8°C
- **Need for Special Timing for Collection (where applicable):** NA
- **Specimen Stability Information:** NA
- **Special Instruction:** NA
- **Rejection Criteria:** Wrong Tube, QNS, Wrong Identification
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Immunoassy	Daily	Next Day	11:00 am	12:00 am	10:00 am