



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital
Clinical Laboratories

Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Clotting Time	CT	Coagulation	Whole Blood

- **Preparation of the Patient:** NA
- **Type of Collection Container:** Plain glass tubes (2)
- **Types and Amounts of Preservatives or Anticoagulants:** NA
- **Quantity of Specimen to be collected:** 1 ml in each tube
 - **Minimum Quantity:** 01 ml in each tube
- **Transportation Temperature:** NA
- **Need for Special Timing for Collection (where applicable):** Sample is not collected at collection point.
- **Specimen Stability Information:** NA
- **Special Instruction:** Not Required
- **Rejection Criteria:** NA
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Manual (Lee & White method)	Daily	Daily	12:01am	Sample is not collected at collection point	Sample is not collected at collection point