



آغا خان یونیورسٹی ہسپتال

The Aga Khan University Hospital  
Clinical Laboratories

### Specimen Collection and Transportation Information

Test Name	Mnemonic	Section	Specimen Type
Urine Phosphorus	URPHOS	Clinical Chemistry	Urine

- **Preparation of the Patient:** Not required
- **Type of Collection Container:** 24 hour Urine Container
- **Types and Amounts of Preservatives or Anticoagulants:** Boric Acid
- **Quantity of Specimen to be Collected:** 24 Hour Urine
  - **Minimum Quantity:** 5cc
- **Transportation Temperature:** 18-25°C (For 8 Hours) 2-8°C after 8 hours to 24 hours
- **Need for Special Timing for Collection (where applicable):** Not applicable
- **Specimen Stability Information:** Stable for 24 hrs. at 2-8 °C and for longer storage at -20°C
- **Special Instruction:** Not required
- **Rejection Criteria:** Wrong identification
- **Need for Appropriate Clinical Data, When Indicated (Patient history):** Not Required

Methodology	Reporting Scheme	Day Performed	Cut Off Time Main Lab	Cut Off Time Karachi Points	Cut Off Time Out of Karachi Points
Phosphomolybdate/UV End Point	Next working day Cutoff 12:00 pm	Daily	12:00 pm	12:00 pm	12:00 pm