آغت خان يونيور سطى به ب يتال

The Aga Khan University Hospital Clinical Laboratories

Specimen Collection and Transportation Information

| Test Name | Mnemonic | Section | Specimen Type |
|---------------|----------|------------|---------------|
| Methemoglobin | METH | Hematology | Whole Blood |

- **Preparation of the Patient:** No special preparation are needed
- Type of Collection Container: Purple top
- Types and Amounts of Preservatives or Anticoagulants: EDTA 7.2 mg
- Quantity of Specimen to be Collected: 4ml
 - o Minimum Quantity: 2-3 ml
- Transportation Temperature: NA
- Need for Special Timing for Collection (where applicable): NA
- Specimen Stability Information: Sample should be as fresh as possible
- **Special Instruction:** Test should not be taken at collection points within and outside Karachi as fresh sample is required for the procedure. Mix the sample well

• Rejection Criteria:

- a) Unlabeled or mislabeled
- b) Identification mismatched between specimen and requisition slip.
- c) Hemolysis, clotted
- d) Improper collection container i.e., leakage or broken sample tube or container
- e) Improper preservative, or improper anticoagulant or medium in container.
- f) Incomplete test requisition slip (i.e., required test information not present)
- g) Specimen volume inadequate (QNS) for analysis
- h) Specimen is from unacceptable source or wrong sample type
- i) Wrong or inappropriate source given for testing
- j) Excessive delay in specimen transport or improperly transported (i.e., not on required temperature)
- k) Specimen contaminated with biological hazardous material
- Need for Appropriate Clinical Data, When Indicated (Patient history): NA

pg. 1

Approved by Service Line Chief and Director Dr. Farooq Ghani dated 1st April, 2016 Controlled Document Original/Copy#: 01/2016

آغت خان يونيور سطى به پتال The Aga Khan University Hospital **Clinical Laboratories**

Specimen Collection and Transportation Information

| Test Name | Mnemonic | Section | Specimen Type |
|---------------|----------|------------|---------------|
| Methemoglobin | METH | Hematology | Whole Blood |

| Methodology | Reporting Scheme | Day Performed | Cut Off Time Main Lab | Cut Off Time Karachi Points | Cut Off Time Out of Karachi Points |
|--------------------------------|---------------------|------------------|--------------------------|---|---------------------------------------|
| Manual (spectrophotometric) | Daily, next day | Daily | 00:01 | Not accepted on collection points | Not accepted on collection points |

pg. 2